

## Using QCM to introduce two or more mutations in one reaction

Sebastian Grünberg, Hahn lab

October 2014

(e.g. changing a His-tag position from N-term. to C-term, multiple nucleotide mutations at two different positions on a plasmid, double deletions, etc.)

This method is based on a 2-step PCR. In the first step, the region between the two regions to be mutated is amplified. In the second step, the whole plasmid gets amplified, using the PCR product from the first step as primers.

This method works surprisingly good, however, careful primer design is important for the QCM to work. For addition of tags (e.g. a 6-His tag), having at least 18 nt upstream and ~ 25 nt downstream of the tag sequence are usually sufficient (if their GC content is about normal). Primers with multiple nucleotide substitutions have to vary in length, acc. to the amount of changed nucleotides. I usually try to have ~ 30 matching nucleotides in the primers.

I use KOD extreme as DNA polymerase for these reactions. However, several other DNAPs (like Phusion, NEB Q5, etc.) work just as well.

One 50 µl rxn contains:

25 µl 2x KOD reaction buffer  
30 ng of plasmid  
250 nM fw primer  
250 nM rv primer  
10 µl 2 mM dNTPs  
1 µl KOD extreme  
H<sub>2</sub>O to 50 µl

For a ~ 680 bp gene on a ~ 5 kb (total) plasmid, I use this program:

1. 95 C 2min
2. 98 C 1 sec
3. 55 C 15 sec
4. 68 C 1 min  
repeat steps 2-4 17x
5. 98 C 10 sec
6. 68 C 14 min  
repeat steps 5-6 17x
7. 68 C 6 min
8. 4 C ∞

EtOH precipitate PCR rxn, resuspend in 30 µl H<sub>2</sub>O  
DpnI digest as in the standard QCM protocol

EtOH precipitate and resuspend in 25  $\mu$ l H<sub>2</sub>O  
Transform 1  $\mu$ l into 25  $\mu$ l of DH10B cells