
Cloning and mutagenesis (TPCR)

Transfer a DNA fragment from plasmid 1 to plasmid 2

1) Design primers:

- Fwd = \sim 20-30 bp overlap with pasmid 2(A) + \sim 20-30 bp overlap with 5'-fragment (from plasmid 1)
- Rev = \sim 20-30 overlap with 3'-fragment (from plasmid 1) + \sim 20-30 bp overlap with pasmid 2(B) and reverse complement.
- If possible, avoid A or T at the extremes of the overlapping regions of the primers. Each of these regions (2 per primer) should have a Tm \sim 60-70 ($Tm = (\sum GC)*4 + (\sum AT)*2$)



2) Prepare reactions (25ul reaction is best):

| Reagent | Cc | Volume (ul) |
|------------|----------|-------------|
| plasmid 1 | 10 ng/ul | 1 |
| plasmid 2 | 10 ng/ul | 1 |
| primer Fwd | 1 uM | 1 |
| primer Rev | 1 uM | 1 |
| dNTPs | 10 mM | 1 |
| Buffer | 5x | 10 |
| H_2O | | 34.3 |
| Phusion | 2U/ul | 0.65 |
| TOTAL | | ~50 |
| | | |

3) run PCR

| temp | time | |
|------|-------------|-----|
| 95 | 60" | |
| 95 | 30" | x35 |
| 60 | 60" | x35 |
| 72 | 4' * | x35 |
| 72 | 7' | |
| 7 | ∞ | |

* → 5 min or longer times for larger plasmids.

4) DpnI digest (100% efficiency of DpnI in Phusion buffer): 10-25ul reaction + DpnI for 90min to O.N. @37°C

5) Transform into E.coli