

## Quickchange Mutagenesis with KOD Extreme

October 2011 Hahn Lab  
B. Knutson/S. Hahn

Quickchange Mutagenesis with KOD extreme DNA Polymerase  
This is much faster than using Pfu Polymerase and typically gives a high frequency of the desired mutation. Works for point mutations, small insertions, and deletions. For large insertions, we use PCR products rather than oligonucleotides (see below)

### **Each reaction contains:**

25 microliters 2X KOD extreme buffer  
30 ng plasmid  
0.2 micrograms ea oligonucleotide (both top and bottom strands ~50-80 bases long depending on mutation)  
10 microliters 2 mM dNTP mix  
H<sub>2</sub>O to 49 microliters total  
1 microliter KOD extreme (Novagen)

1. 98 deg 2 min
2. 98 deg 1 sec
3. 55 deg 10 sec
4. 68 deg 1 min/KB

repeat steps 2-4 for a total of 18 times.

5. Hold at 4 deg

*(The mutagenesis reaction can probably be scaled down to 25 microliters volume but we haven't yet tried this)*

Ethanol precipitate the mutagenesis reaction and resuspend DNA in 30 microliters H<sub>2</sub>O.

Digest the mutagenesis reaction with Dpn1 to remove the starting plasmid DNA:

10 microliters of PCR reaction from above  
16 microliters H<sub>2</sub>O  
3 microliters NEB buffer 4  
1 microliter Dpn1 enzyme

Digest 3 hrs at 37 deg (can do for 2 hr if pressed for time)

Transform 1 microliter of Dpn1 digested DNA to E. coli and select for transformants.

Typically get >50% mutagenesis efficiency.

### **Generating large insertions with Quickchange:**

For large insertions, we use PCR products instead of top and bottom oligonucleotides. Design the oligos for PCR amplification with ~35 nucleotides at the 5' ends that are identical to the template sequence flanking the insertion point.

Do the PCR amplification in a 50 microliter reaction. For the mutagenesis to work, the PCR reaction must work very well (a very prominent band on a gel when ~5% of the total reaction is run on an agarose gel). A minor band or two is OK - no need to purify the PCR product.

Ethanol precipitate and resuspend the amplified DNA in 30 microliters H<sub>2</sub>O.

Perform quickchange with KOD extreme as above, except use 13 microliters of the PCR product that was resuspended in H<sub>2</sub>O instead of oligonucleotides.

Ethanol precipitate, and digest with Dpn1.

Typically get >50% insertions.