

## READY TO GO DNA SEQUENCING GUIDELINES

**Tube/Plate requirements:** Samples should be submitted numerically in strip tubes with strip lids. We recommend strip tubes from Fisher cat#14-230-210.

Plates must be loaded A-H with no spaces. V-bottom plates only, half or non-skirted.

**Recommended DNA Quantities:** DNA should be in dH<sub>2</sub>O or Tris buffer. Please ensure DNA templates are free of salts, proteins, organics, and florescent tags. PCR products should be free of extra primers and dNTP's.

Template Type	DNA Quantity/Reaction (ng)
<b>PCR products</b>	
100 to 200 bp	0.5 to 3
200 to 500 bp	1 to 10
500 to 1000 bp	2 to 20
1000 to 2000 bp	5 to 40
>2000 bp	10 to 50
<b>Plasmids</b>	
<6 kb	50 to 300
6-10 kb	300 to 500
>10 kb	500 to 1,000
<b>Other types of template</b>	
single-stranded	10 to 50
Cosmid or BAC DNA	200 to 1,000
Bacterial genomic DNA	1,000 to 3,000

### Preparing Reactions:

Reaction size	BigDye mix	Primer (pmoles)	Template	DMSO** (5-10% final volume)	H <sub>2</sub> O to final volume	Number of cycles
1/2 reaction in 10ul	4ul	3-5	See above	0.5-1ul	10ul	25
1/2 reaction in 20ul	4ul BD + 4ul 5X buffer	3-5	See above	1-2ul	20ul	25
1/4 reaction in 10ul	2ul	3-5	see above	0.5-1ul	10ul	25
1/8 <sup>†</sup> reaction in 10ul	1ul BD + 1.5ul 5X buffer*	3-5	see above	0.5-1ul	10ul	25
1/8 <sup>†</sup> reaction in 5ul	1.0ul BD + 0.5ul 5X buffer*	3-5	see above	0.25-0.5ul	5ul	25

Note: Applied Biosystems full reaction is 8ul.

\*5x sequencing buffer is an Applied Biosystems product and is available in DE-302 at no additional charge

\*\*DMSO is an optional reagent and will not hurt the reactions. It is useful for GC-rich sequences.

†1/8 reactions should be tried only after success at 1/4 reaction

**Thermal Cycling Conditions:** These thermal cycling conditions were optimized using GeneAmp PCR system 9700, Applied Biosystems 9800 Fast Thermal Cycler, the Veriti 96-Well Thermal Cycler, and the Veriti FAST 96-Well Thermal Cycler. If you choose other thermal cyclers, you may need to adjust the conditions.

**For Double-stranded DNA, Single-stranded DNA, and PCR Products:**

Stage	Description	Temp (°C)	Time
1	Denaturation	96	1 min
2	Amplification: 25 cycles	96	10 sec
		50	5 sec
		60	4 min
3	Hold	4	Indefinite hold

**For BAC DNA:**

Stage	Description	Temp (°C)	Time
1	Denaturation	95	5 min
2	Amplification: 50 cycles*	95	30 sec
		50-55**	10 sec
		60	4 min
3	Hold	4	Indefinite hold

\*Some laboratories have found increasing number of cycles gives better results.

\*\*Set the annealing temperature according to the template.

**Additional information:**

-For templates with secondary structure we recommend using dGTP BigDye. Aliquots are available in the DE-302 freezer.

-For small templates and reading close to the primer we suggest decreasing both BigDye and template.

-In-House primers are M13 Forward, M13 Reverse, T7, and SP6. These are in the DE-302 refrigerator at 3uM concentration and are no additional charge.

**Note from Hahn lab: For sequencing double stranded plasmids, use the conditions: 1/8 reactions with 10 microliters total volume. This works well for nearly everything and saves on Big dye reagent – Steve Hahn.**