

Sporulation and Dissection of Yeast a/alpha Diploids

Linda Warfield

Last Modified Tue, Nov 6, 2001

These conditions work well for the Hieter yeast strains

1. Isolate diploid colonies. Inoculate 2 ml YPD + Ade with a single, large diploid colony and grow O/N at 30 degrees C to stationary phase.
2. Put 200 ul of stationary phase cells into 5 ml sterile, distilled H₂O. Spin down 2 min. in clinical centrifuge. Wash again with 5 ml H₂O.
3. Resuspend in 2 ml 0.5% (50 mM) KOAc (pH 7.0) + 0.5X nutrients for auxotrophic markers (usually Ade, Ura, Trp, Leu, His, Lys, Met).
4. Incubate on tube roller at rm. temp. for 3-7 days.
5. Check for sporulation under microscope. Spin down tetrads in a small glass tube and wash 3 times with 5 ml sterile H₂O. Resuspend in 2 ml H₂O.
6. Spin down 180 ul cells. Resuspend in 90 ul Zymolyase, 0.5 mg/ml in 1M sorbitol. Incubate at 30 degrees for 5-10 minutes.
7. Slowly add 0.3 ml H₂O on ice to stop the reaction.
8. Plate 30 ul cells onto YPD + Ade, Trp, Ura and dissect.