

Standard Operating Procedure 2.5: Preparing Cryopreserved Cell Suspension for Transplant

1. Thaw frozen vial of cryopreserved tumor cells in warm water bath (~37°)
2. Immediately after thawing is complete, add 15ml of serum-free media to dilute the DMSO that the cells are frozen in
3. Determine cell number and cell viability
 - a. Pipet 1 ml of cell/media mixture into a ViCell-specific cuvette
 - b. Place cuvette in the ViCell counting tray and press “Log in Sample”
 - c. Enter a sample ID, set cell type to “default” and dilution factor to “1”
 - d. Press “Start Queue” to begin counting
4. With the information from the ViCell, calculate the resuspension volume using this formula:

Resuspension volume (in microliters) =

$$\frac{\text{Viable Cells/ml (Vicell reading)} \times \text{total vol. of cell suspension (ml)} \times 1000 \text{ ul/ml}}{50 \times 10^6 \text{ cells/ml}}$$

5. Spin the remaining 15 ml cell suspension at 1000 rpm for 5 minutes
6. After spin has finished, discard all supernatant media using a vacuum pipet, leaving just the cell pellet in the conical
7. Resuspend the pellet in the volume of serum free media determined in step 8, which will be at a final concentration of 5×10^6 cells per ml (100,000 cells/2ul)
8. Transfer cell suspension to a 1-2 ml tube and place on ice
9. Fill out the logsheet in the tissue culture room

*NOTE: if you are also making a cell culture dish use this same cell suspension (assuming there are enough cells) and refer to SOP3 for plate and media preparation.