# **Chromatin IP (Isw2)**

7/01 Toshi last update: 06/15

### **Reagents**

Fix soln: 11% formaldehyde, 0.1 M NaCl, 1 mM EDTA, 50 mM Hepes-KOH pH 7.6. <u>Freshly prepared</u>. Do not store in glass bottles.

2.5 M glycine.

TBS: 20 mM Tris pH 7.6, 150 mM NaCl.

Breaking buffer: 100 mM Tris pH 8.0, 20% glycerol, 1 mM PMSF (added fresh).

Acid washed glass beads 425-600  $\mu$  (Sigma G-8772).

FA buffer: 50 mM Hepes-KOH pH 7.6, 150mM NaCl, 1 mM EDTA, 1% Triton X-100, 0.1% sodium deoxycholate. NaCl conc must be determined for each combination of tagged factor and epitope tag. Add protease inhibitors (100x) just before use. [100X cocktail: 0.1 M PMSF, 0.2 mM pepstatin A, 60  $\mu$ M leupeptine, 0.2 mg / ml chymostatin, 0.2 M benzamidine in 100% EtOH. Store at -20 °C].

FA-HS buffer: same as FA buffer but contains 500 mM NaCl.

- RIPA buffer: 10 mM Tris pH 8.0, 0.25 M LiCl, 0.5% NP40, 0.5% sodium deoxycholate, 1 mM EDTA.
- 3XFLAG peptide (Met-Asp-Tyr-Lys-Asp-His-Asp-Gly-Asp-Tyr-Lys-Asp-His-Asplle-Asp-Tyr-Lys-Asp-Asp-Asp-Asp-Lys. 5 mg / ml stock in Buffer H 0.1.) Crude peptide is sufficient.

Dynabeads Protein G (Dynal, cat# 100.04).

- FLAG M2 antibody (Sigma).
- 0.1 M sodium Phosphate buffer pH7.0
- 0.1 M sodium Phosphate buffer pH7.0 + 0.01% Twwn 20

2x Stop buffer: 20 mM Tris pH8.0, 100 mM NaCl, 20 mM EDTA, 1% SDS.

### **Chromatin Prep**

- 1. On the day before chromatin prep, inoculate yeast so that the culture will be  $OD_{600}=0.5$  at the time you want to start chromatin prep. I routinely prep chromatin from 50-300 ml culture.
- 2. At OD<sub>600</sub>=0.5, add 1/10 volume of Fix soln to culture, and shake slowly @ room temp for 20 min. <u>Optimum time for fixation must be determined for</u> <u>each epitope tag and tagged protein</u>. We have done 5 or 20 min fix for Isw2 and saw the same results.
- 3. Add 18 ml of 2.5 M glycine to every 100 ml of culture to terminate fixation. Shake culture for 5 min at room temp. If you are taking time course, you can keep the cells on ice for a while (up to a couple of hors) at this stage.
- 4. Pellet cells (5,000 rpm 5 min in JA-10), and wash twice with ice cold TBS (10 ml ea) in 14 ml tubes. You can freez cell pellet in liquid N2 and store in -80 °C at this stage.
- 5. Re-suspend cells in breaking buffer (500  $\mu$ l buffer for 300 ml culture) and transfer to flat bottom 2 ml microfuge tubes. Jairo uses 300  $\mu$ l breaking buffer for cells from 50 ml culture.
- 6. Add glass beads (500  $\mu$ l for 300 ml culture, ~100  $\mu$ l for small scale) and bead beat 2-5 minutes. <u>It's very important not to over bead beat</u>. <u>Check cells</u> <u>under microscope</u>, and stop when 70-80 % of cells are broken. Optimum time widely varies among different strains.
- Puncture one hole each in the bottom then top of the tube by a needle, and put tubes on Falcon 2059 (14 ml snap cap) or equivalent tubes. Centrifuge @ 1,000 for 2 min to recover cell lysates.
- 8. Add 1 ml FA buffer to lysates, mix, and transfer to 2 ml microfuge tubes, and spin at 14,000 rpm for 1 min at 4 °C and discard sup.
- 9. Suspend ppt in 1 ml FA buffer by gentle pipetting. Repeat step 8 twice to partially purify chromatin. suspend chromatin in 1 ml FA buffer.
- 10 Sonicate up in Biorupter. I normally do 15 min sonication twice at 50% cycle. For some prep, 30 minutes sonication is too much. If so, reduce it to 25 minutes total.
- Centrifuge in microfuge tubes at 14,000 rpm for 30 min twice (turn tubes in between) at 4 °C, and transfer the sup into a new eppi tube. The sup is the chromatin sample. Go to IP reaction, and/or freeze chromatin in liquid N<sub>2</sub>. Store at -80 °C.

# Beads Prep

1. Take 20  $\mu$ I of Pro G beads/sample (IP reactions) into a siliconized eppi tube.

- 2. Concentrate the beads on a concentrator (MPC) and remove sup. <u>To avoid</u> <u>clamping, do not leave beads on MCP for more than a couple of minutes</u>.
- 3. Add 0.5 ml of 0.1 M sodium phosphate pH7.0 (or PBST).
- 4. Rotate at room temp for 2~3 min. Concentrate on MPC and remove sup.
- 5. Repeat step 4 twice (total 3 washes).
- 6. Take up beads in 20~100  $\mu$ I 0.1 M sodium phosphate pH7.0 (or PBST) and suspend beads. Add anti-FLAG Ab (normally ~2  $\mu$ I/IP: need to be optimized for each protein). If multiple IP is to be done, steps 6-8 should be done in one tube.
- 7. Constantly mix at room temp for 30 minures.
- 8. Concentrate on MPC and remove sup. Add 0.5 ml 0.1 M sodium phosphate pH7.0 + 0.01% Tween 20 (or PBST), suspend beads and rotate 2~3 minutes. Concentrate on MPC. Repeat beads wash x1 in 0.5 ml FA biffer, and suspend beads at 20  $\mu$ I/IP. Keep beads on ice until use.

#### **IP** reaction

- 1. Thaw chromatin in room temp water, centrifuge at the top speed for 5 minutes and take the amount needed. I usually take 240  $\mu$ l (200  $\mu$ l for IP and 20  $\mu$ l each for DNA and western analyses).
- 2. Save 20  $\mu$ l each of chromatin for DNA and western blotting ("input" sample). Add 20  $\mu$ l of 2 x Stop buffer or 2 x SDS/PAGE sample buffer to DNA and western samples, respectively.
- 3. Add 200  $\mu$ l of chromatin from step 2 to 20  $\mu$ l of Protein G Dynabeads as prepared above in siliconized microfuge tubes. Rotate at room temp for 90 min for binding.
- 5. Concentrate beads on MPC. Save 20  $\mu$ l of sup for western ("unbound" fraction) and add 20  $\mu$ l 2 x SDS/PAGE sample buffer.
- 6. Add 1 ml FA buffer, suspend beads by flicking and rotate 5 min at room temp. Concentrate on MPC and remove the buffer.
- 7. Do washing (step 6) for total 3 times with FA buffer, 2 times with FA-HS buffer and 1 times with RIPA buffer.
- 8. Elute by adding 45  $\mu$ l of x2 Stop buffer and incubating at 65 °C for 15 minutes with occasional vortexing.
- 12. Save 5  $\mu$ l of eluate for western from each eluate ("bound" fraction). Add 20  $\mu$ l 2 x SDS/PAGE sample buffer and 15  $\mu$ l H2O.
- 13. For western blotting, incubate samples for 30 min on 95 °C heat block. Determine how much of Isw1/2 was precipitated by comparing signals for "input", "unbound" and "bound".

- 14. For all DNA samples, add 1  $\mu$ l of 20 mg / ml glycogen per tube. Add 20  $\mu$ l of 2x Stop buffer and 20  $\mu$ l H2O (total 80  $\mu$ l). Incubate @ 75 °C for at least 6 hrs, then add 2  $\mu$ l of 10 mg / ml proteinase K. Incubate at 50-55 °C for several hrs to overnight. Phenol-chloroform extract twice and EtOH ppt DNA (20  $\mu$ l 10 M NH4OAC and 250  $\mu$ l EtOH).
- 14\*: Alternatively, you can prepare DNA by Qiagen PCR Purification Kit. Elute in 50  $\mu$ l.
- \*Ultra-rapid DNA prep method [ref: Nelson et al., Nature Protocols <u>1</u>:179-185, 2006).
- Required reagent: 10 % (wt/vol) Chelex 100 (Bio-Rad, cat # 142-1253) in distilled water.
- 1. Input DNA (20  $\mu$ l at the step 3 above): add 60  $\mu$ l EtOH. Cfg at max for 15 minutes and wash ppt twice with 80% EtOH. Suspend DNA in 100  $\mu$ l 10 % (wt/vol) Chelex 100 and dissolve DNA.
- 2. Bound DNA: after the last wash (RIPA at the step 7 above), add 100  $\mu$ l 10 % Chelex 100 directly to the beads. Vortex.
- 3. Boil both input and bound DNA 30 minutes on 95 °C heat block. Take 10  $\mu$ l for western blotting ("bound" sample).
- 4. Add 2  $\mu$ I 10 mg/ml PK to each tube and incubate 55 °C 30 min to overnight. I have done 30 minutes, a few hours and overnight, and all worked fine.
- 5. Boil 10 minutes on 95 °C heat block. Vortex, cfg and recover sup.
- 6. Add 100  $\mu$ l water to beads, vortex, cfg and take the sup. Combine samples from steps 5 and 6, add 1  $\mu$ l glycogen and EtOH ppt. Rinse ppt twice with 80 % EtOH and dissolve them in 50  $\mu$ l TE.

## PCR

- 1. Take up DNA in 50  $\mu$ l of TE. Take 3  $\mu$ l of " input " samples, treat with RNase A and run on 1.0-1.3 % agarose gel to check size of DNA. You should see smears centered around 300-500 bp.
- 2. For control, I usually use 5-fold serial dilution of " input " sample. Try 10, 50, and 250-fold dilution. For "bound", use undiluted sample.
- 3. Typical PCR reaction (20  $\mu$ l).

2 μl DNA sample 1X buffer 1.5 mM MgCl<sub>2</sub> 20 pmole each primers 200 μM each 4dNTP

95 °C 2 min--[95 °C 20 sec - 59 °C 40 sec - 72 °C 30 sec] x 26-28 cycles--72 °C 5 min (RoboCycler).

- 4. Run on 2 % agarose in 0.5 x TBE.
- \*:You may need to try different number of cycles for different primers/samples. For Isw1/2, I adjust amount of DNA so that 26 cycles would be sufficient to see signals. Higher numbers of cycles are not recommended because reactions tend to go out of linear range.
- \*\*: For radioactive PCR for quantitation (strongly recommended), I add 0.5  $\mu$ I [α-<sup>32</sup> P] dCTP (3,000 Ci/mmol, 10  $\mu$ Ci/ $\mu$ I) per 100  $\mu$ I PCR reaction. Run DNA on 2% agarose in 0.5 x TBE.
- \*\*\*: For CHIP on DNA microaray, I follow Rich Young lab's protocol (http://staffa.wi.mit.edu/cgibin/young\_public/navframe.cgi?s=12&f=ChIPArray) after step 1.